

Film Array Standard Standard for Polariscope

It is available for color sample of retardation, and also, for standard piece to calibrate a Polariscope.

FAS300

- ❖ Convenient to quantify retardation in a sample
- ❖ Available to watch color gradation by retardation on Sensitive tint color method
- ❖ Including with the certification of test.



Specification

Measurement range : 0-200nm
Accuracy : ± 5 nm
Material : Coatable Retarder
Size : 45x210x10mm
Weight : 100g

The retardation values on above photo are reference. Each product has own values.

Orihara Industrial Co., Ltd.
5-47-15, Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013, Japan
TEL +81-3-3985-9531 FAX +81-3-3985-9532
<http://www.orihara-ss.co.jp> e-mail orihara@orihara-ss.co.jp